TRETTS




Gebrauchsanleitung
User's Manual

Komponenten:

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Components:

1. Main frame
2. Front boom
3. Gear hub
4. Under-seat steering
5. Seat
6. Twist shifter
7. Brake lever
8. Seat adjustment
9. Frame number / Bike coding
10. Battery headlight
11. Battery rear-light
12. Parking support
13. Drawbar
14. Front wheel bracket

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Scope of delivery:

- Yellow green RAL 6018
- Frame, front boom, fork aluminium
- Seat height 13 - 17 in. (33 - 43cm)
- Seat recessed in the middle for safe positioning of the trunk
- Seat can be adjusted to the bottom bracket position without tools
- Lower handlebars can be positioned in width and inclination
- Tires Kenda 58 - 406mm with Reflective strips and puncture protection
- Spoke protectors
- Crankset Alu shortened to 115 mm to reduce the pedal amplitude 32Z
- 8 speed Nexus premium gear hub
- SKS Mudguards black matt
- Lighting system with rechargeable battery: Spanninga Axendo headlight, Spanninga Elips taillight
- Five-point harness
- Children's pedals with toe clips and straps heel rubber for safe positioning of the feet
- non telescopic safety flag

In this basic configuration, the TRET5 is already equipped for the disabled. Further information on p. 24

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Introductory notes

Contraindications

The TRETS is not suitable for humans with visual impairment.

Further information

If you have any problems or questions about your therapeutic bike or trike, please contact your authorized dealer or HASE BIKES for assistance.

Any serious incident* that occurs in connection with the product must be reported to the manufacturer and the responsible competent authorities of the member state in which the user and/or patient resides.

* A "serious incident" is any incident that had, could have had, or could have one of the following consequences, either directly or indirectly: a) the death of a patient, user, or other person, b) the serious temporary or permanent deterioration of the health condition of a patient, user, or other person, c) a serious risk to public health.

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


Adjusting your TRETS

The proper adjustment of the TRETS is important for riding comfort and performance. The TRETS can be adjusted for use by riders from approx. 1.00 to 1.60 m (3'2" to 5'2") in height. Thanks to the adjustable seat, the TRETS can be adjusted without tools by 15 cm. In addition to the length adjustment, it is also possible to individually adjust the seat height to your own requirements.

For larger adjustments an adaptation of the frame and the chain is necessary. Please refer to the following pages how to adjust the seat height and the frame length.

Checking the seat position

Sit down on the TRETS seat and place your feet on the pedals with the balls of your feet on the pedal axles. If one pedal is at the farthest point, the leg should be almost straight and in a comfortable position. If the pedals are too close (e.g. it is not possible to extend the leg) or too far (e.g. the rider cannot easily reach both pedals), then the length of the TRETS must be adjusted.

Adjusting the seat position

1. Loosen the three levers (1) on the front seat mount.
2. Move the seat to the desired position and tighten the previously released levers.
3. Check the seat position. If any changes are necessary, repeat the steps.

Note:

If you want to adjust the length by more than approx. 15 cm, you have to adjust the frame length.

Adjusting the seat inclination

1. Loosen the two quick releases (2) on the rear seat attachment.
2. Move the seat to the desired position and tighten the two quick releases.
3. Secure the position additionally by clipping in the seat support clips in the highest possible position (3) on both sides.

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


Adjusting the seat height

Adjusting the seat height at the mounting plate

By means of the mounting plate (1) the seat height can be increased by a total of 5cm in two steps. Loosen the clamping levers and turn the mounting plate to the desired position.

Adjusting the seat height at the seat

The seat height can be increased by an additional 8 cm on the seat. Unscrew the clamping lever on the seat (2) and change to the lower hole.

Adjusting the frame length

If the adjustment range of the seat is extended or the desired seat angle is not reached due to an extreme size adjustment, the frame length can be set individually.

Adjusting the frame length to the rider additionally optimizes the center of gravity of the trike and thereby improves the cornering stability, the turning circle and the riding comfort.

Note:

A frame length change is only possible in combination with a shortening or extension of the chain. In addition to a chain tool, you will also need the appropriate chain links for an extension. The TRETS comes with an 9-speed chain.

1. First loosen the M6 screw on the clamp (3) of the handlebar.
2. Unscrew the four bolts (4) at the rear frame (see page 5).
3. Push the front boom slightly backwards to take tension off the chain. Open the chain at the master link and remove the chain.
4. Move the front boom to the desired position.
5. Shorten or lengthen the chain according to the setting. Mount the chain.

Note:

Pull the front boom back slightly to close the master link. The chain tension is then adjusted by a slight push of the front boom.

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

- Tighten the four Allen bolts (2) on the rear triangle with a 5mm Allen key.

Caution!

It is absolutely necessary to tighten these bolts with a torque of 18 Nm (13 ft-lbs). Riding with loose bolts poses a high risk of accident and damage to the frame.

- Adjust the steering. Align front wheel and steering with each other and tighten the M6 screw.

Note:

Any restrictions in handlebar movement could lead to accidents with possible injuries.

Caution!

It is absolutely necessary to tighten these bolt with a torque of 9 Nm (6.6 ft-lbs). Riding with loose bolts poses a high risk of accident.

- Adjust the handlebars and the seat.
- Check the position of the brake and shift cables. Check that the handlebar has full freedom of movement by moving it as far as it will go to the left and to the right.

Adjusting the steering on the TRETS

You can adjust the handlebar width and tilt to your shoulder width and arm length.

- Unscrew the screws (1) on each side of the handlebar with a 5mm Allen key.
- Adjust the inclination of the handlebars.
- Tighten the screws.

Caution!

The handlebar ends may only be pulled out to the marking. Make sure that the handlebars cannot bump against the wheels or fenders. Major risk of accident! The steering could lock in the handlebar's maximum left or right position.

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Riding your TRETS

Correct tire pressure

The TRETS should only be ridden with the correct tire pressure. The tire pressure range is specified on the tire sidewalls. The two rear wheels are equipped with spoke guards. Open the zipper to access the valve.

Note:

Tire pressure affects the handling and safety of the trike. Therefore, always ride with the tire pressure indicated on the side of the tire.

Shifting

Shifting gears allows you to adjust pedaling strength to the terrain and desired speed. The force required per pedal revolution can be decreased by selecting the appropriate gear. In this way, it is much easier to ride up a steep hill in a lower gear than in a high gear. In a lower gear, however, you are slower and have to pedal faster. High gears are best suited for riding at higher speeds on flat terrain or riding downhill. Select a gear that allows you to pedal at more than 70 pedal revolutions per minute. This cadence is optimal for producing sustained power and preventing excessive strain on the knee joints.

The TRETS is equipped with an 8-speed Shimano Nexus hub gear as standard. With a hub gear, you can shift while pedalling lightly or even when the pedals are not in motion.

Note:

Please note the information provided in the enclosed manual from the gear system manufacturer.

Note:

Before riding the TRETS on public roads, be sure to familiarize yourself with gear shifting in an open area away from traffic. Shifting gears takes practice, which could distract you from navigating road traffic.

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

Shifting gears with the twist shifter

Turn the twist shifter (1) on the right handlebar clockwise to shift to a higher gear for faster riding or counterclockwise to shift to a lower gear for slower riding.

Braking

In its standard version, the TRETS is equipped with mechanical disc brakes. The brake levers are mounted on the left and right grips of the handlebar.

The left wheel brake is controlled with the left brake lever, and the right wheel brake with the right brake lever.

Note:

The disc brakes of the bike are very powerful and have good modulation. However, it is still possible for the wheels to become locked. Before riding the TRETS on public roads, be sure to familiarize yourself with the braking force by trying a few braking trials in an open area away from traffic.

Adjusting the brake levers to hand size

To position the brake lever closer to the handlebar grip, turn the adjusting bolt (2) clockwise using a 2mm Allen key. To position the brake lever further from the handlebar grip, turn the bolt counterclockwise.

Parking brake

The TRETS is equipped with a parking brake to prevent it from rolling when parked. There is an additional lever on the brake lever which fixes the brake.

1. Press the brake lever and fix it with the parking brake (3).
2. The brake is now engaged.
3. Press the release lever (4). The built-in spring releases the brake.

Note:

Be sure to release the parking brake before you start.

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


Adjusting the drawbar

The drawbar must be adjusted to the height of the coupling on the towing bicycle before starting the journey.

1. Select the hole on the drawbar connection according to the wheel size of the rear wheel.
2. Insert the M6 screw through the hole of the Deichsel and the drawbar connection and tighten the screw with 9 Nm.

Fit the drawbar

1. Guide the hitch (2) over the front locking pins.
2. Turn the drawbar so that the spring-loaded safety pins lock into the rear hitch of the drawbar.
3. Hitch the trailer to the towing bicycle and remove the front wheel from the trailer.
4. You can fix the front wheel in the holder provided behind the seat.

Caution!

Before starting the journey, check the screw connection of the drawbar connection and the drawbar for a firm fit. Risk of accident!

Caution!

The maximum payload in trailer operation is 45kg.

Note:

The enclosed Weber E-Coupling is only approved for bicycles and pedelecs (max. 25km/h), not for S-Pedelecs or E-Bikes above 25km/h!

Remove the drawbar

1. Loosen the spring-loaded locking pins by pulling the holding wire (3) backwards.
2. Rotate the drawbar and take it off.
3. You can clip the drawbar into the holders provided behind the seat.

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Harness System

For all persons who are at risk of slipping or tilting from the seat while driving or during braking manoeuvres, the use of the 5-point harness included in the delivery is recommended.

Mount the lap belts to the seat as shown in the picture and the shoulder belts to the diagonal braces of the rear seat connection. The fifth strap can be attached to both the stem and the lower seat tube.

Caution!

When riding, the buckles of the Harness System must be securely fastened at all times. Otherwise, they could get caught in moving components.

Note:

Loose ends of straps can be cut off using a hot cutter. Consult your authorized dealer.

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Pedal with toe Clips and straps


This is the simplest and most universal solution for securing the foot in the pedal. The straps can be individually adjusted to the size of the foot.

Adjustment and Use

The toes are inserted into the toe clips (1) and the elastic band (2) is pulled around the heel. Adjust the length of the toe-clip strap and the elastic band to fit the foot. The length of the elastic band can be adjusted by changing the position of the knot in the band.

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Flag

A non-telescopic flag is included in the delivery.

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Parking the TRETS upright

The TRETS can be parked upright with the front wheel raised. In this position, it requires less space for parking than a standard bicycle. To bring the TRETS into the upright parking position, simply lift the front wheel.

Note:

Due to the adjustable seat and the frame length, there is a danger of over-tipping when the settings are extreme. Make sure you have a safe place. Such accidents could cause damage to your TRETS or other nearby objects, such as bicycles, cars, etc.

If the TRETS has a Steerable Push Bar (available as accessory), the trike cannot be parked upright on its rear wheels.

Transporting the TRETS on a car rack

The TRETS can be transported with a standard roof or rear-mounted bicycle carrier. The rack must be equipped for three bicycles and adjusted so that each TRETS wheel can be anchored on a separate wheel tray. Secure all wheels and the frame using load straps. Check the mounting to be sure that the TRETS cannot come loose during the trip. Remove all parts that could fall off in transport.

Transporting the TRETS

The TRETS can be adjusted to its most compact state for transportation by car, train, or other vehicle. The front boom can be fully inserted and the seat can be removed easily.

1. Loosen the quick release clamp of the rear seat attachment and the M8 screw of the front seat attachment. Remove the seat.
2. Release the clamp on the handlebar and tilt the handlebar ends backwards.
3. Slide the seat plate to the bend in the stem.
4. Loosen the front M8 screw (1) of the chain deflection unit.
5. Release the handlebar clamp and adjust the frame to the shortest length. A description can be found in chapter "Adjusting the frame length" on page 4.

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
7. Remove the chain from the front chainring. Make sure that the chain does not damage the frame.
8. Tilt the seat supports mounted on the main frame forwards.
9. Loosen the quick release clamp on the fork and remove the front wheel.
10. Now the TRETs is ready for transport.

The assembly is carried out in reverse order.

Separate the front boom

If the TRETs is still too long, you can separate the front boom with a few further steps

11. Loosen the steering rod clamp. Disconnect the steering rod with two 10mm wrenches by holding the front steering rod (1) and turning the rear steering rod counter-clockwise. Secure the two sections of the steering rod on the front boom and main frame respectively.
12. Loosen the four M6 counter sunk bolts (2) on the rear section of the front boom and separate the two frame parts.





If you are planning to take the trike on an airplane or have it shipped by a forwarding company, we recommend packaging it in a cardboard box for protection. Our TRETs shipping box is very well suited for this purpose.

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Using the light system

Taken from the manufacturer's instructions:

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
Starrer Antrieb
Fixed hub

Nexus Hub – switchable to fixed gear

This hub can be used in two different modes:

A: 8 speeds with freewheel.

B: Fixed hub (helps eliminate the dead spot of the pedal stroke).

i Note: For the TRETs, the only purpose of the 27T sprocket on the hub flange is for holding the ends of the pins (2).

Activating/deactivating the fixed hub.

Switching procedure: Always rotate the shifting ring (1) forwards in the direction of the arrows until the next stop position is reached.

- If the ends of the pins (2) are protruding from the hub flange, the 8 speeds and freewheel are available. (Fig. A)

- If the ends of the pins (2) are not protruding from the hub flange, the fixed hub is activated. (Fig. B)

i **Note:** The cranks must be rotated slightly after switching modes so that the pins can lock into place.

! Important! For safety reasons, the Switchable Hub can only be switched between modes when the trike is parked

! Important! For safety reasons, the fixed hub must be deactivated when the trike is being used in trailer mode. Therefore, check whether the pins (2) on the hub flange protrude approx. 10mm (Fig. A).

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Maintenance intervals

Note:

Please have all inspection and maintenance work performed by your authorized Hase retailer.

The first inspection should be done after the initial 300 km (200 miles). All additional inspections should be performed in 3000 km (2000 mile) intervals or at yearly intervals (see also page 27, "Inspection pass").

Component	Action	Before every ride	Once a month	Once a year
Tires	Check tire pressure Check tread wear and sidewalls	X X		
Brakes	Check pad thickness Check alignment to the brake disc Test braking effectiveness	X	X X	
Brake cables	Check for damage	X		
Chain	Check for wear, oil		X	
Crank bracket	Check bolts			X
Wheels	Check for lateral/radial true Check spoke tension		X X	
Handlebar	Inspect			X
Steering rod	Check for play			X
Headset	Check for play			X
Hubs	Check for play			X
Pedals	Check for axle play			X
Derailleur cables	Check for damage		X	
Quick releases	Check for tightness	X		
Bolts and nuts	Check for tightness		X	
Bottom bracket	Check for play			X

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
Maximaler Luftdruck
Maximum pressure

Checking wheels for lateral/radial true

The wheels have been carefully trued prior to shipment. During the break-in period, stress on the wheels may cause the spokes to settle into the rims so that the spoke tension is no longer uniform and the wheel untrue. Have your wheels re-trued by your authorized retailer at the first inspection after approx. 300 km (200 miles). Subsequent tensioning is seldom necessary. Check the spoke tension and lateral/radial true of the wheels regularly.

You can check the spoke tension by squeezing together two adjacent spokes, working all the way around the rim. If you notice major differences in spoke tension, the wheel should be re-trued. You can check whether a rim is true by attaching a wire to the fork or rear frame and positioning it so that it almost touches the rim. Lift the TRETS and rotate the wheel. If portions of the rim move away from or come in contact with the wire, the wheel must be taken to a qualified bicycle mechanic for re-truing.

Note:

The rear wheels do not have to be removed for the truing process.

Checking tires

Check the tires for damage and tread wear. Tires with cuts, cracks, or worn tread should be replaced. Check tire pressure before each ride and always ride with the pressure specified on the tire sidewalls.

Note:

Never ride with excessive tire pressure (higher than specified on the sidewall). The tire could burst. Risk of accident and injury!

Changing tires or repairing a flat

The front wheel of the TRETS is mounted with a quick release. Simply release it to remove the front wheel.

You do not have to remove the back wheels for changing a tire or repairing a flat. The tires can be removed and remounted using a standard tire lever with the wheels attached.

When removing and remounting tires, take care not to damage the fenders and check that the rim tape is positioned correctly on the rim.

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Checking the disc brakes

Please note the information provided in the enclosed manual from the brake manufacturer. Have a qualified bicycle mechanic check the pad thickness regularly. If you notice any unusual braking sounds or inadequate braking strength, consult your bicycle mechanic immediately.

Checking the chain

The trike's chain is longer than that of a standard upright bicycle. A long chain is generally subject to less wear. Therefore, it does not have to be changed as often as shorter chains.

You can check the condition of the chain using a chain wear gauge from a bike shop.

We recommend that you have a qualified bicycle mechanic check the condition of the chain and, if necessary, replace the chain, and possibly also the chainring and sprocket.

Cleaning and oiling the chain

In order to ensure that the chain does not wear out before the end of its service life, it must be cleaned and lubricated on a regular basis.

1. Pour a small amount of suitable chain cleaner onto a clean, lint-free cotton cloth, and use it to wipe down the chain. When doing so, slowly turn the crank backwards (i.e. opposite the riding direction). Repeat this procedure with clean patches of the cotton cloth until the chain is clean. Allow the cleaning agent to evaporate for approx. 1 hour before applying a lubricant.

Note:

Never use solvents or pressure cleaners. These could wash the lubricants out of the chain bearings and shorten the life of the chain. Never allow chain lubricant to come in contact with the brake discs or tires.

This could significantly reduce braking performance or, in the case of lubricant on tires, wheel traction.

2. Apply chain oil, grease, or wax evenly to the chain links while slowly rotating the crank.

Note:

Rain can wash off traces of lubricant, releasing it into the environment. We recommend biodegradable lubricants.

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3. Rotate a few turns to allow the lubricant to spread.
4. Wait several minutes to allow the lubricant to penetrate. Then wipe off all excess lubricant using a clean, dry cloth.

Note:

If the chain has reached its wear limit, a simultaneous replacement of the cassette is recommended.

Checking the chain guard tube

The purpose of the chain tube (1) is to protect your clothing from chain lubricant. Over time, friction from the moving chain will cause the tube to wear. Check the ends of the chain tube for any signs of wear. If the ends are completely worn, have your bicycle mechanic replace the chain tube.

Note:

In order to replace the chain tube, the chain must be opened at the master link.

Cleaning and disinfecting your TRETS

Be sure to clean your TRETS on a regular basis. Dirt and grime, as well as salt from winter roads, can cause damage. To clean your TRETS, use a sponge and soapy water (ideally with a biodegradable household cleaning product). Stubborn marks on metal surfaces with a matt finish can be removed best by using a sponge with a wheel-rim cleaner from an automotive-supply shop.

You can also clean your TRETS with plain water. If you use a hose, be sure to spray the water gently. Never use pressure washers. The high pressure could force water into the bearings, damaging them. Wipe all surfaces with a moist, lint-free cloth.

Apply a suitable spray wax or a similar corrosion inhibitor to all surfaces of the TRETS except for the brake pads, brake discs, grips, brake/gear levers, chain, seat cover, and tires. Then clean and lubricate the chain. The handlebar grips and seat can be cleaned and/or disinfected using standard cleaning agents and disinfectants. Be sure to use a disinfectant that is suitable for plastic materials.

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Tightening torques

Attention!

Use a torque wrench for all screws. Please check the sa-fety-relevant screws and observe the following table with torques.

Bolt	Tightening torque
M8 lens head bolts at the front boom	13.3 ± 0.7 ft-lbs
M8 Allen bolt in the rear hub	14.7 ± 0.7 ft-lbs
M8 lens head bolts at the fork	13.3 ± 0.7 ft-lbs
M6 counter sunk bolt at the frame seperation	5.9 ± 0.7 ft-lbs
M8 counter nut at the steering rod	11 ± 0.7 ft-lbs

The following table lists the standard values for Allen bolts of 8.8 quality (Grade 5)

Bolt	Tightening torque
M 5	4 ft-lbs
M 6	6.6 ft-lbs
M 8	15.5 ft-lbs

Accepting delivery

If you notice the packaging has been damaged, please notify the shipping company directly.

Unpacking and checking contents

Check the contents of the package with your order.

In addition to the front boom, rear frame, and the seat, the package should also include:

Pedals, decal sets, drawbar, TRETS User's Manual as well as manuals from other manufacturers for various components such as gear system and brakes.

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
Assembling the TRETs

1. Mount the seat. Slide the rear seat attachment over the rear seat tubes (1).
2. Fasten the seat in the desired height (see page 4) with the clamping lever to the seat plate on the front boom.
3. Adjust the position and angle of the seat. Secure the position additionally with the seat support clips (see page 4).
4. Adjust the handlebar tilt (see page 5).
5. Correct the air pressure in the tires.
6. Check all screws for tight fit.
7. Mount any accessories that you have ordered.

Note:

Note that the enclosed chain is designed for an average body size. Due to the slidable seat, the TRETs can be adjusted in this setting to a body size of approx. 1,25-1,40 m (4'1"-4'7").

Shorten or lengthen the chain for a different setting.

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


Removing the wheels

1. Loosen the self-locking M6 nuts (1) of the protective cover with a spanner. Pull the protective cover to the side.
2. Loosen the brake (2) on the adapter.
3. Hold the wheel firmly and turn the M8 screw (3) counter clockwise using a 6 mm Allen key.
4. Remove the wheels together with the brake disc.

Installing the wheels

1. Turn the axle so that the feather key (4) points upwards. Position the groove (5) in the hub to the feather key and slide the wheel together with the brake disc onto the axle.
2. Tighten the wheel with the washer and the M8 screw (3).
3. Fasten the brake disc cover

Attention!

It is absolutely necessary to tighten these bolts with a torque of 21 Nm (17 ft-lbs). Use medium strength loctite. Riding with loose bolts poses a high risk of of accident and damage to the frame.

Dismantle hub gears

1. Hold the wheel firmly and turn the M8 screw (3) counter clockwise using a 6 mm Allen key.
2. Loosen the grub screws at the adjusting collar (6) and the shaft attachment (7) on the right axle.
3. Relax the flange bearings by loosening the four self-locking M6 nuts (8).
4. Slide out the right axle until the nut is free under the shaft attachment. Use a draw hammer if the shaft cannot be pulled out by hand.

Note:

In order to avoid damage to the drive shaft, the draw hammer threaded rod must be screwed into the thread as far as it will go. Alternatively, the brake calliper can be removed and the wheel used as a pulling aid.

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



5. Repeat steps 1 through 4 on the left. On the left side, also loosen the clamp of the hub connection (9).
6. Loosen the screws (10) on the rear frame and push the front boom slightly backwards to remove tension from the chain. Remove the chain from the chainring.
7. The hub can now be removed and the sprocket replaced.

Mounting the drive shaft

1. Push the shaft connection onto the ball bearing and the hub connection onto the hub.
2. Push the drive shaft (6) into its end position. Note the orientation of the feather key (7) to the wheel hub and the setscrews on the shaft to the hub driver (4) or shaft connection.

Attention!

It is absolutely essential that the M8 bolts for securing the hub be tightened with a torque of 21 Nm (15.5 ft-lbs). In addition, medium-strength threadlocker must be used. If the bolts are tightened with insufficient torque, they could come loose, and the drivetrain system could fail.

3. Press the set collars (1) firmly against the outer flange bearing and tighten the set screws.
4. Tighten the M8 screw (3) in the wheel hub.
5. Tighten the four M6 lock nuts (8) on the flange bearings.
6. Tighten the grub screws on the hub driver and shaft connection. Tighten the M6 screws of the clamp to 9Nm.
7. Mount the chain and push the stem slightly to tighten the chain.
8. Tighten the four hexagon socket screws (10) on the rear frame.

Attention!

It is absolutely necessary to tighten these bolts with a torque of 18 Nm (13 ft-lbs). Riding with loose bolts poses a high risk of of accident and damage to the frame.

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Safety information

Intended use

The TRETs is a recumbent trike intended for use on paved roads and bicycle paths, as well as on roads or trails surfaced with sand, gravel, or similar materials (e.g. forest paths, dirt roads). Your safety on roads, trails, and off-road terrain will depend on your riding speed. The higher your speed, the greater your risk! The use of the TRETs for competition or any form of racing is strictly forbidden.

Period of use

With proper use, taking into account the instructions for use and depending on the frequency of use, the expected useful life is up to eight years.

Before your first ride...

Your authorized retailer has assembled your new TRETs and adjusted it to your size requirements. Please refer to the section "Adjusting your TRETs" on page 3 if you need to make adjustments for another rider. Due to its unique construction, there are significant differences in handling in comparison with a standard bicycle. Before riding the TRETs on public roads, it is important to familiarize yourself with the handling in a quiet, open area, away from traffic.

The maintenance measures required before each ride are to be performed by the user (individuals with no physical, mental, or visual impairments).

Observing traffic regulations

The Trets comes with a lighting system and all safety features required by the German road traffic regulations (StVZO) and may only be used in road traffic with this equipment. Be sure to study the applicable regulations and consult your authorized retailer before making any changes to the lighting system.

Traffic regulations are different in every country. The Trets and its equipment must be in compliance with the regulations of the country in which it is used.

Avoiding muscle strain

In a recumbent position, the rider uses different muscle groups than with a standard bicycle. Riders with little or no experience on recumbent bikes/trikes should begin with short, easy trips in order to avoid excessive muscle strain.

Tipping and overturning

Recumbent tricycles have a low center of gravity. This allows for cornering at higher speeds. The centrifugal force felt when turning at high speeds, however, can cause the TRETs to skid or even tip over, which could cause injury to the rider or others. Avoid making hectic steering movements or braking in curves. This increases the TRETs's tendency to tip. Be careful when stopping on inclines. The TRETs could roll backwards and, if brakes are applied abruptly, even overturn. Always wear a helmet for your own safety.

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Braking

The disc brakes of the TRETS are very powerful and have good modulation. However, it is still possible for the wheels to become locked. Before your first trip in public traffic, be sure to take the time to familiarize yourself with the braking strength by trying a few braking trials in an open area away from road vehicles.

Suitable clothing and helmet

Do not wear loose clothing that could get caught in the wheels or gears. Beware of scarves and clothing strings/cords that hang down. Always wear a helmet in order to avoid head injury in the case of an accident.

Risk of injury from rotating parts

If a patient seated in a TRETS that is being pulled as a trailer allows his/her hands or fingers to hang down beside or below the seat, there is a risk that they could get caught in the disc brakes and/or wheels. To prevent the associated injuries, it is essential to ensure that patients keep their hands/fingers on the handlebar grips or above the seat at all times!

Pedals

In a recumbent position, it is especially easy for the feet to slip from the pedals. Always wear shoes with good tread. Maintain sufficient pressure on the pedals. Use toe clips and straps (slip protection) or clipless pedals.

Being visible for others

The TRETS is lower than an upright bicycle. This may make it more difficult for other road vehicles to see you. Always ride with a flag. It attracts the attention of other road users.

Riding at night

Always use lights when riding at night or in conditions of low visibility, such as dusk and dawn.

Technical condition

The TRETS should only be ridden in perfect technical condition. Check for damage before each trip. Take the TRETS to your bike shop for regular maintenance. The recommended maintenance intervals can be found in the section "Maintenance intervals" on page 13.

Replacement parts and accessories

Use only the replacement parts and accessories recommended by HASE BIKES for the TRETS. If other parts are used, the warranty will lose its validity. HASE BIKES is not liable for any resulting damage. Consult your bike shop for more information.

Loading the TRETS

Be sure not to exceed the maximum permissible combined weight for rider and luggage: 100 kg (220 pounds). All luggage should be secured using the mount behind the seat or a lowrider rack (accessory). Always ensure that the load is distributed evenly.

Risk of flammability

The upholstery material of the seat cover is not resistant to ignition sources as specified in EN 1021-2. Therefore, be sure to keep the seat away from open flames, heat sources, etc.

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Description

What is the TRET5?

The TRET5 is a recumbent trike and an approved medical device as per the European Directive (2017/745) MDR (Medical Device Directive). The TRET5 is driven by leg power via the pedals, which are specially designed to prevent the feet from slipping out when riding. It is available with an 8-speed internal-gear hub. On models with the optional differential, both rear wheels are driven, and on standard models, only one rear wheel is driven. The padded seat has been specially designed and constructed with flexible materials to cushion the ride and offer maximum comfort. The TRET5 has a reliable brake system with two independently operated disc brakes.

Indications for use

The TRET5 is suitable for riders with neuromuscular disorders, such as cerebral palsy or muscular dystrophy, who – owing to special physical conditions – are unable to ride standard bicycles or vehicles, even with training wheels. Other indications include coordination disorders, balance disorders, and restricted movement associated with spasticity, paraplegia, dysmelia, amputations, etc.

Contraindications

The TRET5 is not suitable for riders with visual impairment.

Adaptability to individual needs

The TRET5 can be equipped to accommodate riders with various physical impairments. A current list of compatible accessories including instructions can be found on our homepage www.hasebikes.com

The following are just a few examples of available accessories: steerable push bar, crank shorteners, special pedals with or without calf support, pedal pendulums

Multiple users

TRET5 that have been used in accordance with their intended use can be reused by a new rider after the TRET5 has undergone an inspection by an authorized retailer.

Protection against environmental influences

The TRET5 is made from special materials that protect it against corrosion: powder-coated aluminum frame parts, plastics, and stainless steel.

Cleaning and disinfection

The handlebar grips and seat can be cleaned and disinfected using household cleaning agents and disinfectants. For more information, see the section “Cleaning and disinfecting your TRET5” on page 17.

Applicable standards and regulations

The TRET5 has been designed and manufactured in compliance with the following regulations and standards:
European Directive (2017/745) MDR (Medical Device Directive)

Frame number and bike coding

The frame number can be found on the end of the main frame (see illustration on the inside cover of this manual). This is also the only place where bike coding can be engraved. Bike coding helps authorities locate the vehicle owner in the case of theft.

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Conditions of warranty and liability

The statutory warranty applies. The dealer must sell you the HASE BIKE in a safe and road-worthy condition. In addition, the company HASE BIKES offers a warranty of three years (from the date of purchase) against frame breakage in accordance with the following conditions. By registering your bike at <https://hasebikes.com/209-1-Bicycle-passport.html> you can extend this warranty to 5 years.

- HASE BIKES will only be held liable for defects resulting from faulty frame manufacturing and assembly. For assembly parts, the relevant manufacturer is liable.
- This additional warranty only applies to the original owner and to usage in accordance with the user's manual.
- Within the warranty period of 3 years, HASE BIKES will repair, or if necessary replace, any broken frame part free of charge. Any defective parts replaced by HASE BIKES become the property of HASE BIKES.
- Any further claims beyond the scope of this warranty are excluded. HASE BIKES will not cover the costs of assembly or disassembly (e.g., by the retailer) or shipping. We will not accept postal deliveries with insufficient postage.
- The maximum permissible load (rider and luggage) of 100kg (220 lbs) specified for the TRETS must not be exceeded.
- If the TRETS is used in competition or subjected to abnormally excessive stress, the statutory and additional warranty becomes void.
- Damage caused by accidents is not covered by the warranty.
- The warranty becomes void if the TRETS has been fitted with incompatible accessories or if accessories are installed incorrectly. Therefore, always have components replaced and installed by an authorized retailer.
- The warranty becomes void if improper modifications or alterations (e.g., grinding, drilling, bending, etc.) are made to the frame.
- It is not possible to extend the warranty or commence a new warranty period.
- Warranty claims can only be handled by the authorized HASE BIKES retailer from whom the TRETS was purchased or by one of our distribution partners.

Disposal instructions

At the end of its useful life, the product should be disposed of properly, in accordance with the local waste disposal regulations (i.e. not as household waste).

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Have your vehicle inspected by your bicycle dealer at regular intervals. During the high season (March through September), it is best to schedule any required repairs or inspections in advance in order to avoid unnecessarily long downtimes and waiting periods.

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Checklist:

- Check pedals for tightness
- Check headset
- Check handlebar for tightness
- Check front and rear brakes
- Check front and rear wheel bearings
- Check bottom bracket
- Check cranks for tightness
- Check seat for tightness
- Check chain tension and wear
- Check tires and air pressure
- Check rims and spokes
- Check all nuts and bolts for tightness
- Check gear system
- Check accessories

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